

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

#### www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB08868 - Goat Anti-TFAP2D (aa70-82) Antibody

Size: 100µg specific antibody in 200µl



#### **Target Protein**

**Principal Names:** TFAP2D, transcription factor AP-2 delta (activating enhancer binding protein 2 delta), TFAP2BL1, AP-2 like, OTTHUMP00000016593, activating enhancer binding protein 2 beta-like 1, transcription factor AP-2 beta (activating enhancer binding protein 2 beta)-like 1, transcription factor AP-2 beta (activating enhancer-binding protein 2 beta)-like 1, transcription factor AP-2 beta-like 1

Official Symbol: TFAP2D

Accession Number(s): NP\_758438.2

Human GeneID(s): 83741

Non-Human GeneID(s): 226896 (mouse)

### **Immunogen**

Peptide with sequence C-HHQSFHYEFQHSH, from the internal region of the protein sequence according to NP\_758438.2.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

#### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:2000.

## **Species Reactivity**

Tested:

Expected from sequence similarity: Human, Mouse