



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB08868 - Goat Anti-TFAP2D (aa70-82) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** TFAP2D, transcription factor AP-2 delta (activating enhancer binding protein 2 delta), TFAP2BL1, AP-2 like, OTTHUMP00000016593, activating enhancer binding protein 2 beta-like 1, transcription factor AP-2 beta (activating enhancer binding protein 2 beta)-like 1, transcription factor AP-2 beta (activating enhancer-binding protein 2 beta)-like 1, transcription factor AP-2 beta-like 1

**Official Symbol:** TFAP2D

**Accession Number(s):** NP\_758438.2

**Human GeneID(s):** [83741](#)

**Non-Human GeneID(s):** 226896 (mouse)

### Immunogen

Peptide with sequence C-HHQS FHYEFQ HSH, from the internal region of the protein sequence according to NP\_758438.2.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:2000.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse