

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05536 - Goat Anti-Thymosin-like 8 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TMSB15A, TMSNB, TMSL8, thymosin, beta, identified in neuroblastoma cells, Thymosin-like 8, thymosin beta 15a, TMSL8, Tb15, TbNB

Official Symbol: TMSB15A

Accession Number(s): NP_068832.1

Human GeneID(s): [11013](#)

Immunogen

Peptide with sequence SDKPDLSEVEKFDR-C, from the N Terminus of the protein sequence according to NP_068832.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: No signal obtained yet but low background observed in 293 and A549 lysates at up to 2µg/ml (please note that our SDS PAGE/western blotting protocol is not optimized for detection of proteins smaller than 15kDa) .

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat