

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05536 - Goat Anti-Thymosin-like 8 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** TMSB15A, TMSNB, TMSL8, thymosin, beta, identified in neuroblastoma cells, Thymosin-like 8, thymosin beta 15a, TMSL8, Tb15, TbNB

**Official Symbol:** TMSB15A

**Accession Number(s):** NP\_068832.1

**Human GeneID(s):** [11013](#)

### Immunogen

Peptide with sequence SDKPDLSEVEKFDR-C, from the N Terminus of the protein sequence according to NP\_068832.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:64000.

**Western blot:** No signal obtained yet but low background observed in 293 and A549 lysates at up to 2µg/ml (please note that our SDS PAGE/western blotting protocol is not optimized for detection of proteins smaller than 15kDa) .

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Rat