

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10559 - Goat Anti-TIAM2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TIAM2, T-cell lymphoma invasion and metastasis 2, FLJ41865, STEF, TIAM-2, OTTHUMP00000017472, SIF and TIAM1-like exchange factor, T-lymphoma invasion and metastasis-inducing protein 2

Official Symbol: TIAM2

Accession Number(s): NP_036586.2; NP_001010927.1

Human GeneID(s): [26230](#)

Non-Human GeneID(s): 24001 (mouse)

Important Comments: This antibody is expected to recognize all reported isoforms (NP_036586.2; NP_001010927.1).

Immunogen

Peptide with sequence C-QIRHQSLDSQSEN, from the internal region of the protein sequence according to NP_036586.2; NP_001010927.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave an approx 70kDa band in NIH3T3 lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 192kDa according to Mouse NP_036008.2. The 70kDa band was successfully blocked by incubation with the immunizing peptide.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Dog