



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05339 - Goat Anti-TIM3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TIM3, T cell immunoglobulin mucin-3, KIM-3, FLJ14428, kidney injury molecule-3, hepatitis A virus cellular receptor 2

Official Symbol: HAVCR2

Accession Number(s): NP_116171.3

Human GeneID(s): [84868](#)

Immunogen

Peptide with sequence QQPSQLGCRFAMP, from the C Terminus of the protein sequence according to NP_116171.3.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 30-35kDa band seen in Human Brain lysate. Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human