

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05339 - Goat Anti-TIM3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TIM3, T cell immunoglobulin mucin-3, KIM-3, FLJ14428, kidney injury

molecule-3, hepatitis A virus cellular receptor 2

Official Symbol: HAVCR2

Accession Number(s): NP_116171.3

Human GenelD(s): 84868

Immunogen

Peptide with sequence QQPSQPLGCRFAMP, from the C Terminus of the protein sequence according to NP_116171.3.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 30-35kDa band seen in Human Brain lysate. Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human