



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08298 - Goat Anti-TMPRSS4 (aa245-57) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TMPRSS3, transmembrane protease, serine 4, MT-SP2, membrane-type serine protease 2, transmembrane serine protease 3, type II membrane serine protease

Official Symbol: TMPRSS4

Accession Number(s): NP_063947.1; NP_001077416.1

Human GeneID(s): [56649](#)

Important Comments: This antibody is expected to recognise both reported isoforms (NP_063947.1 and NP_001077416.1).

Immunogen

Peptide with sequence HCFRKHTDVFNWK, from the internal region of the protein sequence according to NP_063947.1; NP_001077416.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Approx. 48kDa band observed in lysates of cell lines HEK293 and A549 (calculated MW of 48.2kDa according to NP_063947.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB08298 (1 μ g/ml) staining of HEK293 lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.