

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08298 - Goat Anti-TMPRSS4 (aa245-57) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TMPRSS3, transmembrane protease, serine 4, MT-SP2,

membrane-type serine protease 2, transmembrane serine protease 3, type II membrane

serine protease

Official Symbol: TMPRSS4

Accession Number(s): NP_063947.1; NP_001077416.1

Human GeneID(s): 56649

Important Comments: This antibody is expected to recognise both reported isoforms

(NP_063947.1 and NP_001077416.1).

Immunogen

Peptide with sequence HCFRKHTDVFNWK, from the internal region of the protein sequence according to NP_063947.1; NP_001077416.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Approx. 48kDa band observed in lysates of cell lines HEK293 and A549 (calculated MW of 48.2kDa according to NP_063947.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB08298 (1 μ g/ml) staining of HEK293 lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.