



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11280 - Goat Anti-TSPYL2 (aa2-13) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CASK-interacting nucleosome assembly protein, CDA1, cell division autoantigen 1, CINAP, CTCL, CTCL tumor antigen se20-4, CTCL-associated antigen se20-4, cutaneous T-cell lymphoma-associated antigen se20-4, cutaneous T-cell lymphoma-associated tumor antigen se20-4, DENTT, differentially expressed nucleolar TGF-beta1 target, differentially-expressed nucleolar TGF-beta1 target protein, HRIHFB2216, NP79, nuclear protein of 79 kDa, SE204, testis-specific protein Y encoded-like 2, testis-specific Y-encoded-like protein 2, TSPX, TSPY-like 2, TSPY-like protein 2, TSPYL2

Official Symbol: TSPYL2

Accession Number(s): NP_071400.1

Human GeneID(s): [64061](#)

Immunogen

Peptide with sequence DRPDEGPPAKTR-C, from the N Terminus of the protein sequence according to NP_071400.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Pig, Cow