



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB08151 - Goat Anti-UBE1L2 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** UBE1L2, ubiquitin-activating enzyme E1-like 2, FLJ10808, FLJ23367, MOP-4, UBA6, UBA6, ubiquitin-activating enzyme E1

**Official Symbol:** UBE1L2

**Accession Number(s):** NP\_060697.4

**Human GeneID(s):** [55236](#)

### Immunogen

Peptide with sequence C-ETTEVRKTKIRN, from the internal region of the protein sequence according to NP\_060697.4.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:4000.

**Western blot:** Western Blot: Preliminary experiments in Human Testis and Liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human