



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05766 - Goat Anti-UBE2F/NEDD8-conjugating enzyme Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** UBE2F, ubiquitin-conjugating enzyme E2F (putative), MGC18120, NCE2, NEDD8 conjugating enzyme, NEDD8-conjugating enzyme

**Official Symbol:** UBE2F

**Accession Number(s):** NP\_542409.1

**Human GeneID(s):** [140739](#)

**Important Comments:** This antibody is expected to recognise isoform 1 (NP\_542409.1) only.

### Immunogen

Peptide with sequence C-RNKVDDYIKRYAR, from the C Terminus of the protein sequence according to NP\_542409.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Preliminary experiments gave no signal but low background in Human Kidney and Liver lysates at up to 1µg/ml.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Rat, Dog, Pig, Cow