

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB09397 - Goat Anti-UHMK1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: UHMK1, U2AF homology motif (UHM) kinase 1, KIS, Kist, KIS protein kinase, OTTHUMP00000029529, kinase interacting stathmin, kinase interacting with

leukemia-associated gene (stathmin)

Official Symbol: UHMK1

Accession Number(s): NP_787062.1

Human GeneID(s): 127933

Non-Human GenelD(s): 16589 (mouse), 246332 (rat)

Immunogen

Peptide with sequence C-AGDSKAAQKLLTGR, from the C-Terminus (near) of the protein sequence according to NP_787062.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:2000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow