

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11706 - Goat Anti-Unc5b (mouse, rat) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: 6330415E02Rik, A630020F16, D10Bwg0792e, netrin receptor UNC5B, protein unc-5 homolog 2, protein unc-5 homolog B, unc5 homolog 2, unc-5 homolog 2, unc-5 homolog B (C. elegans), Unc5b, Unc5h2

Official Symbol: Unc5b

Accession Number(s): NP_084046.2

Human GeneID(s): [219699](#)

Non-Human GeneID(s): 107449 (mouse), 60630 (rat)

Important Comments: This antibody is expected to recognize the cytoplasmic domain of the protein.

Immunogen

Peptide with sequence C-DWIFQLKTQAHQGH, from the internal region of the protein sequence according to NP_084046.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments in Human and rodent Brain, Lung and Spleen lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow