

International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05984 - Goat Anti-USP11 / UHX1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: USP11, UHX1, ubiquitin specific protease 11, deubiquitinating enzyme 11, ubiquitin thiolesterase 11, ubiquitin carboxyl-terminal hydrolase, X-linked, ubiquitin specific peptidase 11, RP4-659F15.2, ubiquitin-specific processing protease 11

Official Symbol: USP11

Accession Number(s): NP_004642.2

Human GeneID(s): [8237](#)

Immunogen

Peptide with sequence CSSPPSSEFMDVN, from the C Terminus of the protein sequence according to NP_004642.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 120kDa band observed in lysates of cell line Jurkat (calculated MW of 110kDa according to NP_004642.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Ideguchi et al, Biochem J. 2002 Oct 1;367(Pt 1):87-95.; PMID: 12084015). Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB05984 staining (1 μ g/ml) of Jurkat lysate (RIPA buffer, 30 μ g total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.