



EB08669 - Goat Anti-USP9X Antibody

Size: 100µg specific antibody in 200µl



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

Target Protein

Principal Names: USP9X, ubiquitin specific peptidase 9, X-linked, DFFRX, FAF, Drosophila fat facets related, X-linked, deubiquitinating enzyme FAF-X, fat facets protein related, X-linked, ubiquitin specific peptidase 9, X-linked (fat facets-like, Drosophila), ubiquitin specific protease 9, X chromosome (fat facets-like Drosophila), ubiquitin specific protease 9, X-linked, ubiquitin specific protease 9, X-linked (fat facets-like, Drosophila), ubiquitin thiolesterase FAF-X, ubiquitin-specific processing protease FAF-X

Official Symbol: USP9X

Accession Number(s): NP_001034679.2

Human GeneID(s): [8239](#)

Important Comments: This antibody is expected to recognize isoform 3 (NP_001034679.2)

Immunogen

Peptide with sequence C-KKATSVQQIEMEES, from the C Terminus (near) of the protein sequence according to NP_001034679.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:16000.

IHC: Immunohistochemistry: In paraffin embedded Human Cortex shows strong cellular body staining of neurons. Recommended concentration, 3-5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB08669 (3.8 μ g/ml) staining of paraffin embedded Human Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.