



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10484 - Goat Anti-VN1R1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: V1RL1, V1R-like 1, V3R-related, VNR19I1, vomeronasal 1 receptor 1, vomeronasal olfactory receptor, (chromosome 19) subtype I, member 1, ZVNH1, ZVNR1, VN1R1

Official Symbol: VN1R1

Accession Number(s): NP_065684.1

Human GeneID(s): [57191](#)

Immunogen

Peptide with sequence C-SSDLNENQHP, from the internal region of the protein sequence according to NP_065684.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Species Reactivity

Tested:

Expected from sequence similarity: Human