



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EBL014-CD68-P2 - Positive Control Lysate for anti- CD68 antibody

Size: 50ul

### Name

Positive Control Lysate for anti- CD68 antibody (EB12699 Batch Specific #P2)

### Amount supplied

50ul (5ug/ul protein in RIPA buffer).

### Test numbers

7

### Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

### Shelf life

3 months.

### Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

### Application tested

Western Blot: 2µg/ml staining of Primary antibody in lysates EBL014-CD68-P2,  
Approx.110kDa band observed (calculated MW of 37.4kDa according to NP\_001242.2).  
Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at  
1:40K



Western Blot: 2 $\mu$ g/ml staining of Primary antibody in lysates EBL014-CD68-P2, Approx.110kDa band observed (calculated MW of 37.4kDa according to NP\_001242.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K