

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EBL014-CD68-P2 - Positive Control Lysate for anti- CD68 antibody

Size: 50ul

Name

Positive Control Lysate for anti- CD68 antibody (EB12699 Batch Specific #P2)

Amount supplied

50ul (5ug/ul protein in RIPA buffer).

Test numbers

7

Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

Shelf life

3 months.

Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

Application tested

Western Blot: 2µg/ml staining of Primary antibody in lysates EBL014-CD68-P2, Approx.110kDa band observed (calculated MW of 37.4kDa according to NP_001242.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

Western Blot: 2µg/ml staining of Primary antibody in lysates EBL014-CD68-P2, Approx.110kDa band observed (calculated MW of 37.4kDa according to NP_001242.2). Primary incubation 1 hour at room temperature.

Secondary antibody is EB2ND-002-A at 1:40K