

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EBL008-CD209-P1 - Positive Control Lysate for anti- DC-SIGN / CD209 antibody

Size: 50ul

Name

Positive Control Lysate for anti- DC-SIGN / CD209 antibody (EB12715 Batch Specific #P1)

Amount supplied

50ul (5ug/ul protein in RIPA buffer).

Test numbers

7

Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

Shelf life

3 months.

Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

Application tested

Western Blot: 2µg/ml staining of Primary antibody in lysates EBL008-CD209-P1, Approx.50kDa band observed (calculated MW of 45.8 kDa according to NP_066978.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K

250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa

Western Blot: 2µg/ml staining of Primary antibody in lysates EBL008-CD209-P1, Approx.50kDa band observed (calculated MW of 45.8 kDa according to NP_066978.1). Primary incubation 1 hour at room temperature.

Secondary antibody is EB2ND-002-A at 1:40K