



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EBL003-GAPDH-S3G2 - Positive Control Lysate for anti- GAPDH antibody

Size: 50ul

Name

Positive Control Lysate for anti-GAPDH antibody (EB06377 Batch Specific #S3G2)

Amount supplied

50ul (5ug/ul protein in RIPA buffer).

Test numbers

7

Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

Shelf life


3 months.

Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

Application tested

Western Blot: 0.1µg/ml staining of Primary antibody in lysates EBL003-GAPDH-S3G2, approx. 37kDa band observed. (Calculated MW of 36.1kDa according to Human NP_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K



Western Blot: 0.1µg/ml staining of Primary antibody in lysates EBL003-GAPDH-S3G2, approx. 37kDa band observed. (Calculated MW of 36.1kDa according to Human NP_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K