



## UK Office

**Everest Biotech Ltd**  
Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## **EBL005-GAPDH-S3 - Positive Control Lysate for anti-GAPDH antibody**

Size: 50ul

### **Name**

Positive Control Lysate for anti-GAPDH antibody (EB06377 Batch Specific #S3G2)

### **Amount supplied**

50ul (5ug/ul protein in RIPA buffer).

### **Test numbers**

7

### **Storage and Handling**

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

### **Shelf life**


3 months.

### **Reconstitution**

Bring vial to 4 °C and add 100ul RIPA buffer

### **Application tested**

Western Blot: 0.5µg/ml staining of Primary antibody in lysates EBL005-GAPDH-S3, approx. 37kDa band observed. (Calculated MW of 36.1kDa according to Human NP\_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K



Western Blot: 0.5 $\mu$ g/ml staining of Primary antibody in lysates EBL005-GAPDH-S3, approx. 37kDa band observed. (Calculated MW of 36.1kDa according to Human NP\_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K