

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EBL012-GPR94-G2 - Positive Control Lysate for anti- GPR94 / TRA1 antibody

Size: 50ul

Name

Positive Control Lysate for anti- GPR94 / TRA1 antibody (EB07290 Batch Specific #G2)

Amount supplied

50ul (5ug/ul protein in RIPA buffer).

Test numbers

7

Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

Shelf life

3 months.

Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

Application tested

Western Blot: 1ug/ml staining of Primary antibody in lysates EBL012-GPR94-G2, Approx. 100kDa band observed (calculated MW of 92.5kDa according to Human NP_003290.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

Western Blot: 1ug/ml staining of Primary antibody in lysates EBL012-GPR94-G2, Approx. 100kDa band observed (calculated MW of 92.5kDa according to Human NP_003290.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K