

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EBL009-GPR94-G2 - Positive Control Lysate for anti- GPR94 / TRA1 antibody

Size: 50ul

#### Name

Positive Control Lysate for anti- GPR94 / TRA1 antibody (EB07290 Batch Specific #G2)

## **Amount supplied**

50ul (5ug/ul protein in RIPA buffer).

#### **Test numbers**

7

# Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

## Shelf life

3 months.

## Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

## **Application tested**

Western Blot: 0.1-0.5µg/ml staining of Primary antibody in lysates EBL009-GPR94-G2, Approx. 100kDa band observed (calculated MW of 92.5kDa according to Human NP\_003290.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K

250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa

Western Blot: 0.1-0.5µg/ml staining of Primary antibody in lysates EBL009-GPR94-G2, Approx. 100kDa band observed (calculated MW of 92.5kDa according to Human NP\_003290.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K