

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EBL018-NUP85-P1 - Positive Control Lysate for anti- NUP85 antibody

Size: 50ul

#### Name

Positive Control Lysate for anti-NUP85 antibody (EB06201 Batch Specific #P1)

## **Amount supplied**

50ul (5ug/ul protein in RIPA buffer).

#### **Test numbers**

7

#### Storage and Handling

Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.

## Shelf life

3 months.

## Reconstitution

Bring vial to 4 °C and add 100ul RIPA buffer

# **Application tested**

Western Blot: 0.1µg/ml staining of Primary antibody in lysate EBL018-NUP85-P1, approx. 35kDa band observed. (Calculated MW of 36.1kDa according to Human NP\_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

Western Blot: 0.1μg/ml staining of Primary antibody in lysate EBL018-NUP85-P1, approx. 35kDa band observed. (Calculated MW of 36.1kDa according to Human NP\_002037.2). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.