
www.everestbiotech.com
Email: customerservice@vectorlabs.com
Telephone: (650) 697-3600

SB-SHIFIX25 COVERSLIPS

SKU: SB-SHIFIX25



DESCRIPTION

Search

Type in product code or name:

SHI-FIX ™ COVERSLIPS

Do you wish immunofluorescence of suspension cells was as easy as with adherent cells? Now it is! Shi-fix™ (manufactured by Shikhar Biotech) cover-slips allow suspension cells to be layered on them or directly grown as a monolayer. Easy, hassle-free protocol and no spin columns needed for fixing suspension cells to slides. Just add your cells to the coverslips, wait for 30 mins, wash unbound cells with PBS and proceed to immunostaining (or continue the culture to obtain desired cell densities).

These cover-slips can also be used for staining suspension cells for microscopy, or for fixing suspension cells for confocal microscopy. What's more, you can also stimulate the cells while they are attached on the cover-slip, should you need to do so. You can work on non-adherent cells as if they were adherent!

Suspension cell immunofluorescence with Shi-fix™ coverslips

SPECIFICATIONS

Unit Size	Pack of 25
Shipping Instructions	Ambient

GALLERY IMAGES

